

The FlexorAide+ It's There When You Need It.

BUY IT ONCE. USE IT A LIFETIME. GUARANTEED. Made In the USA.™/Pat Pending

~ “Few health problems have the potential to affect the horse regardless of age, breed, gender, housing, management, or intended use the way lameness can. All types of horses are susceptible to lameness, and sooner or later, most horses are affected by it.” (Jeffcott 1982, Linder 1993), Equine 98 Census Report.

~ Every year hundreds of thousands of the equine population from recreational to various disciplines will become lame due to injuries, aging, etc. **Flexor tendons** are the most common, with four limbs susceptible to this specific injury, and an **82%** chance of reoccurrence.

The FlexorAide+ is the only product of its design and multipurpose function to be made available for direct purchase. This purpose is to provide **safe, immediate immobilization to prevent further injury.** The FlexorAide+ **promotes continued support** throughout the stages of the regenerative healing cycle with the ongoing supervision of your veterinarian, regarding the appropriate manipulation of the limb and proper setting of the splint during the first 3 to 4 months of regeneration and fibre organization, a crucial point in the rehabilitation and promotion of normal flexion/load bearing.

The FlexorAide+ is simple by design, with a user friendly, single person, 3 step application to front or hind limbs for nonsurgical/post-surgical progressive rehabilitation cases such as hyper extension or rupture of superficial digital/deep digital flexor tendons, suspensory ligaments, condylar fractures, laminitis, proximal suspensory desmitis, minor cannon/splint bone fractures (below knee/hock), abrasions, lacerations, etc.

The FlexorAide+ is 100% RADIOLUCENT, enabling 360-degree radiographs and assessment of the injury during stages of healing, without removal of the splint, providing a safe, practical and easy way to progressively adjust the limb and eliminate the possibility of stress, discomfort, pain, or re-injury.

The FlexorAide+ is made entirely of eco-friendly polycarbonate, one of the strongest materials on the planet, capable of withstanding 9000+ PSI (950 to 1200lb horse, hoof to surface is approx 8000 pounds of pressure per hoof at full gallop), it's padded with carbon black, closed cell, polyethylene foam (enabling radiolucency). It's durable yet provides comfort and the prevention of abrasions. The foam is affixed with **3M 9472LE industrial strength adhesive, formulated to continue adhesion against various weather conditions and aggressive solvents, allowing a continuous wash and use cycle.**

The FlexorAide+ provides an alternative to other devices and procedures regarding flexor tendon injuries, such as, expensive rental fees for “precision instruments” “in need of regular servicing and main component assessment and replacement on regular intervals”.

Or therapeutic shoeing, a work intensive/expensive process where various thickness of wedges are attached to the bottom of the hoof and periodically changed approximately every 4 to 5 weeks (by the farrier) who also keeps the hoof trimmed and balanced. The FlexorAide+ does not involve the hoof to achieve progressive healing, therefore eliminating this arduous process.

The FlexorAide+ is a creation of LMI Direct. A company founded on 23 years of experience in the field of Veterinary medicine. Its founder's hands-on care in the field led him to create safer and more effective tools for the trade, with the intention of providing a less stressful and more practical method for the patient and provider. This would eventually lead to the line of essential products for the veterinary, (and eventually) human medical and dental fields.